



ULTRASHAPE

The Shape of Things to Come

The First Effective Non-Invasive Ultrasonic Solution For Body Contouring



Make a Measurable Difference to More Patients' Lives

- A safe "walk-in walk-out" office procedure
- No recuperation period
- The first effective treatment option when surgery is not desired

Innovative Technology; Selective Effects



- Selective emulsification of adipose tissue by focused ultrasound
- Clearance of the lysed adipocytes via physiological mechanisms
- Skin, connective tissues, blood vessels and nerves remain intact
- Computerized tracking algorithm guarantees homogenous coverage of the treatment area

Where Safety and Effectiveness Finally Meet

- The only treatment that combines the safety of non-invasive body contouring with measurable results
- An average of 2 cm circumference reduction in a single treatment

